

CAS KIT (HCRSNK)





CAS Full KIT (HCRSNKP)






129

OSSTEM KI

CAS KIT Surgical Instruments

CAS Drill

- · Specailized drill designed to penetrate the sinus floor without damaging the schneiderian membrane by forming a concial bone lid
- Four blade body drills well at both high and low speeds and is capable of collecting autogenous bone at low speeds
- Use with stoppers for safe and controlled penetration
- Final drill should be based on the bone quality, regardless of the fixture type (straight or tapered)
- Recommended speed: 400~800rpm (first time: 400rpm)



| $L \setminus D$ | Ø2.8 | Ø3.1 | Ø3.3 | Ø3.6 | Ø3.8 | Ø4.1 |
|-----------------|------------|------------|------------|------------|------------|------------|
| Short | SNDR2813TS | SNDR3113TS | SNDR3313TS | SNDR3613TS | SNDR3813TS | SNDR4113TS |
| Long | SNDR2813TL | SNDR3113TL | SNDR3313TL | SNDR3613TL | SNDR3813TL | SNDR4113TL |

Guide Drill

- Marks the fixture's insertion site
- Side cutting blades trim the extraction socket sidewalls
- Marker 2mm from the tip


Ø2.0/2.7 D SNGD2027TL



Ø2.2 Twist Drill

- Recommend under-drilling by 1mm less than the bone's thickness
- Use with stoppers for safe and controlled drilling
- The tip measures an additional 0.6mm


Ø2.2 \ D SNTD2213TL



Hydraulic Membrane Lifter Set

- Hydraulic pressure is used to separate and lift the sinus membrane
- Securly fits Ø 2.8~ Ø 4.1 CAS drilled osteotomies

\ **D** Ø2.6/6.0 SNMLS



Stopper

- Laser marked numbers indicate the remaining tool's (drill, instruments, etc.) length
- Color-coded by length
- Drill and stopper recommended number of usage is 50 times



Bone Carrier

- Handle for the bone carrier head
- Connect the bone carrier head and tighten at the opposite end
- Connects both heads (SNBCH30 or SNBCH35)





Bone Carrier Head

- Cone shaped with an extended tip that reaches the sinus cavity and prevents bone material from spilling out
- SNBCH30 for Ø 3.1/3.3 CAS drilled osteotomy
- SNBCH35 for Ø 3.6/3.8/4.1 CAS drilled osteotomy
- Fill the reservior with bone material (up to the marker), with the bone condenser shuttle the material in small quantities into the sinus. Repeat the process as necessary.



CAS KIT Surgical Instruments

Bone Condenser

- Safely shuttles bone material through the bone carrier head into the sinus cavity
- SNBCH30 : use Ø 1.1 / SNBCH35 : use Ø 1.4





Hydraulic Membrane Lifter Tube


• Tubing connects to the hydraulic membrane lifter and salin filled syringe

SNMT



Membrane Lifter

- Round shape, no cutting edge and safe membrane lift
- After the CAS drill is used, the membrane was lifted and the lifter diameter was selected according to the CAS drill diameter (head diameter : CAS drill diameter -0.2mm)
- Using CAS stopper for depth adjustment
- Recommended speed: 400~800rpm (for first user: 400rpm)
- Be sure to spray water when using



| L D | Ø2.6 | Ø2.9 | Ø3.1 | Ø3.4 | Ø3.6 | Ø3.9 |
|-------|------------|------------|------------|------------|------------|------------|
| Short | SNML2813TS | SNML3113TS | SNML3313TS | SNML3613TS | SNML3813TS | SNML4113TS |
| Long | SNML2813TL | SNML3113TL | SNML3313TL | SNML3613TL | SNML3813TL | SNML4113TL |

Depth Gauge

• Measures the thickness of the residual bone and checks to see if the sinus is properly separated from the floor


SNDG




Bone Spreader

- A tool that spreads a filled bone by using engine
- Used with stopper
- Recommended speed: 30rpm or less (low speed mode)

Ø2.0 Ø3.0 \ D SNBS2015T SNBS3015T





Y-Connector

• Y-type connecting tool capable of simultaneous water pressure elevation in two drilling holes

SNYCT

