

Digital Library

Distribution of the latest version

2024.11

01 **OSSTEM[®]** x 3shape 
IMPLANT

1. Check the version

- ✓ Ensure that the library remains **current** by updating it

Example
Name)



① Manufacturer ② System ③ Scan Source ④ Target ⑤ Revision



✓ Update Highlights

- Additional features / specification
- Error correction
- Improvement....

2. Check the available libraries

- ✓ Check if the product is being used in local country

Example Name)

OST_TS_SBO_CU_R2

Check !

Check !

③ Scan Source
: Product to be Scanned



④ Target
: Type of prosthesis to be fabricated



Example)







① Manufacturer	② System	③ Source (Product to be Scanned)		④ Target (Type of prosthesis to be fabricated)		⑤ Revision				
Osstem	System	Scan body		Custom Abt.		Revision				
OST/OSSTEM	TS	SB	(Old) Scan body	CU	Customized Abt.	R0	Initial Library File			
		SBO	New Scan Body	Link Abt.		R1	Library Revision 1st			
	LinkSB	Link Scan body (Link Only S.B)	LI_V1	Link (A)	• • •					
			LI_V2	Link (B)						
	SS	MultiSB	Multi Scan body (Multi Only S.B)	Stock Abt.						
		Scan Healing		TR				Transfer Abt.		
		SH/SH(V1)	Scan Healing	Mu-Ti				Multi Ti-Base		
	US	Stock Abt.		MU				Multi Abt.		
		TR	Transfer Abt.	Mu(A)				Multi Angled Abt.		
		CO	ComOcta Abt.	MU_Ti				Multi Abt. + Muti Ti-base		
		COP	ComOcta Plus Abt.	Mu(A)_Ti				Multi Angled Abt. + Multi Ti-base		
	KS	CE	Cement Abt.	CO				ComOcta Abt.		
				COP				ComOcta Plus Abt.		
				CE				Cement Abt.		
OST	TS/SS/US/KS	SB/SH/TR/...		CU/LI/TR/...				R0/R1...		

List of Digital Library : 3 Shape

System	Name of library	Revision	Scan Source		Target		
			Product Name	Code	Product Name	Code	
KS	OST_KS_SBO_CU_R*	4	New Scan Body	KSNSB*	Custom Abt		KSPM*
	OST_KS_SBO_TR_R*	1			Transfer Abt		KSTA*
	OST_KS_SBO_LI_R*	2			Link Abt		KSPL*
	OST_KS_SB_CU_R*	3	Scan Body	KSSB*	Custom Abt		KSPM*
	OST_KS_SB_TR_R*	0			Transfer Abt		KSTA*
	OST_KS_SB_LI_R*	0			Link Abt		KSPL*
	OST_KS_SH_CU_R*	2	Scan Healing	KSSHA*	Custom Abt		KSPM*
	OST_KS_TR_TR_R*	0	Transfer Abt	KSTA*	Stock	Transfer Abt	KSTA*

List of Digital Library : 3 Shape

KS	Scan Source	New Scan body	Scan Body	Scan Healing	Transfer Abt.		
							
		KSNSB*	KSSB*	KSSHA*	KSTA*		
	Target	Custom Abt	Transfer Abt.	Link Abt.			
							
		KSPM*	KSTA*	KSPL*			












List of Digital Library : 3 Shape

System	Name of library	Revision	Scan Source		Target			
			Product Name	Code	Product Name	Code		
TS	OST_TS_SBO_CU_R*	2	New Scan Body	TSNSB*	Custom Abt		TSPM*/*PM*	
	OST_TS_SBO_TR_R*	1			Transfer Abt		GSTA/TSTA*	
	OST_TS_SBO_LI_V1_R*	3			Link	Link Abt (A)		TSPTB*
	OST_TS_SBO_LI_V2_R*	4				Link Abt (B)		TSPL*
	OST_TS_SBO_LI_ASC					Link Abt (B) ASC		TSL*
	OST_TS_SBO_MU_R*	1			Multi	Multi Abt		TSMA/GSMA*
	OST_TS_SBO_Mu_Ti_R*	2				Multi Abt + Ti-base		+ TSMTB*
	OST_TS_SBO_Mu(A)_R*	1			Multi Angled	Multi Angled Abt		GS/TS*MA*
	OST_TS_SBO_Mu(A)_Ti_R*	2				Multi Ang. + Ti-base		+ TSMTB*

List of Digital Library : 3 Shape

System	Name of library	Revision	Scan Source		Target			
			Product Name	Code	Product Name	Code		
TS	OST_TS_SB_CU_R*	4	Scan Body	TSSB*	Custom Abt		TSPM*/*PM*	
	OST_TS_SB_TR_R*	2			Transfer Abt		GSTA/TSTA*	
	OST_TS_SB_LI_V1_R*	3			Link	Link Abt (A)		TSPTB*
	OST_TS_SB_LI_V2_R*	3				Link Abt (B)		TSPL*
	OST_TS_SB_MU_R*	2			Multi	Multi Abt		TSMA/GSMA*
	OST_TS_SB_Mu_Ti_R*	0				Multi Abt + Ti-base		+ TSMTB*
	OST_TS_SB_Mu(A)_R*	0			Multi Angled	Multi Angled Abt		GS/TS*MA*
	OST_TS_SB_Mu(A)_Ti_R*	0				Multi Ang. + Ti-base		+ TSMTB*
	OST_TS_MultiSB_Mu_R*	1	Multi Scan Body	TSMABC	Multi / Multi Angled Abt		TSMA/GSMA* GS/TS*MA*	
	OST_TS_MultiSB_Mu_Ti_R*	1			Multi / Multi Angled Abt + Ti-base		+ TSMTB*	
	OST_TS_SH(V1)_CU_R*	2	Scan Healing	TSSHA*	Custom Abt		TSPM*/*PM*	
	OST_TS_SH(V1)_LI_V2_R*	0		TSSHA*	Link	Link Abt (B)	TSPL* (Shared with TSL*)	
	OST_TS_TR_TR_R*	1	Transfer Abt	TSTA*/GSTA*	Stock	Transfer Abt	TSTA*/GSTA*	
	OST_TS_LinkSB_LI_R*	1	Link Abt Scan body	TSPTB*	Link	Link Abt (A)	TSPTB*	

List of Digital Library : 3 Shape




TS	Scan Source	New Scan body	Scan Body	Scan Healing	Transfer Abt.	Multi Scan Body	Link Scan Body
							
		TSNSB*	TSSB*	TSSHA*	GSTA/TSTA*	TSMSBC	TSPTB*
	Target	Custom Abt	Transfer Abt.	Link Abt (A)	Link Abt (B)	Multi Ti-Base	Multi/ Multi Ang.
							
		TSPM*/*PM*	GSTA/TSTA*	TSPTB*	TSPL* / TSL*	TSMTB*	GS/TSMA* GS/TS*MA*

* TSPL : Link Abutment with EbonyGold screw
* TSL : Link Abutment with Ti-screw

List of Digital Library : 3 Shape

System	Name of library	Revision	Scan Source		Target			
			Product Name	Code	Product Name	Code		
SS	OST_SS_SBO_CU_P_R*	3	New Scan Body	SSNSB*	Custom Abt		SSPM*	
	OST_SS_SBO_CO_R*	1			Stock	ComOcta Abt		SSCA*
	OST_SS_SBO_COP_R*	2				ComOcta Plus Abt		SSCAP*
	OST_SS_SBO_LI_R*	2			Link	Link Abt		SSPL*
	OST_SS_SB_CU_P_R*	2	Scan Body	SSSB*	Custom Abt		SSPM*	
	OST_SS_SB_CO_R*	1			Stock	ComOcta Abt		SSCA*
	OST_SS_SB_COP_R*	1				ComOcta Plus Abt		SSCAP*
	OST_SS_SB_LI_R*	2			Link	Link Abt		SSPL*
	OST_SS_SH_CU_R*	1	Scan Healing	SSSHA*	Custom Abt		SSPM*	
	OST_SS_CO_CO_R*	0	ComOcta Abt	SSCA*	Stock	ComOcta Abt		SSCA*
	OST_SS_COP_COP_R*	1	ComOcta Plus Abt	SSCAP*		ComOcta Plus Abt		SSCAP*






List of Digital Library : 3 Shape

SS	Scan Source	New Scan body	Scan Body	Scan Healing	ComOcta Abt.	ComOta Plus Abt.	
							
	SSNSB*	SSSB*	SSSHA*	SSCA*	SSCAP*		
	Target	Custom Abt	ComOcta Abt.	ComOta Plus Abt.	Link Abt.		
							
	SSPM*	SSCA*	SSCAP*	SSPL*			

List of Digital Library : 3 Shape

System	Name of library	Revision	Scan Source		Target		
			Product Name	Code	Product Name		Code
US	OST_US_SBO_CU_R*	2	New Scan Body	USNSB*	Custom Abt		USPM*
	OST_US_SBO_CE_R*	1			Stock	Cement Abt	CAM/CAR/CAW*
	OST_US_SBO_LI_R*	1			Link	Link Abt	USPL*
	OST_US_SB_CU_R	3	Scan Body	USSB*	Custom Abt		USPM*
	OST_US_SB_CE_R*	0			Stock	Cement Abt	CAM/CAR/CAW*
	OST_US_SB_LI_R*	0			Link	Link Abt	USPL*
	OST_US_SH_CU_R*	1	Scan Healing	USSHA*	Custom Abt		USPM*
	OST_US_CE_CE_R*	0	Cement Abt	CAM/CAR/CAW*	Stock	Cement Abt	CAM/CAR/CAW*

List of Digital Library : 3 Shape

US	Scan Source	New Scan body	Scan Body	Scan Healing	Cement Abt.		
							
		USNSB*	USSB*	USSHA*	CAM/CAR/CAW*		
	Target	Custom Abt	Cement Abt.	Link Abt.			
							
		USPM*	CAM/CAR/CAW*	USPL*			

02 **OSSTEM[®]** *x* **exocad**
IMPLANT

Check the available libraries

- ✓ Check if the product is being used in local country

Example Name)

 **Osstem_TS_SBO_CU** ✓ Check! ✓ Check!

③ Scan Source
: Product to be Scanned



④ Target
: Type of prosthesis to be fabricated



List of Digital Library : Exocad



X exocad

System	Name of library	Scan Source		Target	
		Product Name	Code	Product Name	Code
KS	Osstem_KS_SBO_CU	New Scan Body	KSNSB*	Custom Abt	KSPM*
	Osstem_KS_SBO_TR			Transfer Abt	KSTA*
	Osstem_KS_SBO_LI			Link Abt	KSPL*
	Osstem_KS_SB_CU	Scan Body	KSSB*	Custom Abt	KSPM*
	Osstem_KS_SB_TR			Transfer Abt	KSTA*
	Osstem_KS_SB_LI			Link Abt	KSPL*

List of Digital Library : Exocad

KS	Scan Source	New Scan body	Scan Body	Scan Healing	Transfer Abt.		
							
		KSNSB*	KSSB*	KSSHA*	KSTA*		
	Target	Custom Abt	Transfer Abt.	Link Abt.			
							
		KSPM*	KSTA*	KSPL*			

List of Digital Library : Exocad



X exocad

System	Name of library	Scan Source		Target			
		Product Name	Code	Product Name	Code		
TS	Osstem_TS_SBO_CU	New Scan Body	TSNSB*	Custom Abt		TSPM*/*PM*	
	Osstem_TS_SBO_TR			Transfer Abt		GSTA/TSTA*	
	Osstem_TS_SBO_LI_V1			Link	Link Abt (A)		TSPTB*
	Osstem_TS_SBO_LI_V2				Link Abt (B)		TSPL*
	Osstem_TS_SBO_LI_ASC				Link Abt (B) ASC		TSL*
	Osstem_TS_SBO_MU			Multi	Multi Abt		TSMA/GSMA*
	Osstem_TS_SBO_MuTi				Multi Abt + Ti-base		+ TSMTB*
	Osstem_TS_SBO_MuATi				Multi Ang. + Ti-base		+ TSMTB*

List of Digital Library : Exocad















X exocad

System	Name of library	Scan Source		Target			
		Product Name	Code	Product Name	Code		
TS	Osstem_TS_SB_CU	Scan Body	TSSB*	Custom Abt		TSPM*/*PM*	
	Osstem_TS_SB_TR			Transfer Abt		GSTA/TSTA*	
	Osstem_TS_SB_LI_V1			Link	Link Abt (A)		TSPTB*
	Osstem_TS_SB_LI_V2				Link Abt (B)		TSPL*
	Osstem_TS_SB_MU			Multi	Multi Abt		TSMA/GSMA*
	Osstem_TS_MultiSB_Mu	Multi Scan Body	TMSBSC	Multi / Multi Angled Abt		TSMA/GSMA* GS/TS*MA*	
	Osstem_TS_MultiSB_MuTi			Multi / Multi Angled Abt + Ti-base		+ TSMTB*	
	Osstem_TS_SH_V1_CU	Scan Healing	TSSHA*	Custom Abt		TSPM*/*PM*	
	Osstem_TS_SH_LI		TSSHA*	Link	Link Abt (B)	TSPL* (Shared with TSL*)	
	Osstem_TS_TR_TR	Transfer Abt	TSTA*/GSTA*	Stock	Transfer Abt	TSTA*/GSTA*	
	Osstem_TS_LinkSB_LI_V1	Link Abt Scan body	TSPTB*	Link	Link Abt (A)	TSPTB*	

List of Digital Library : Exocad



X exocad








TS	Scan Source	New Scan body	Scan Body	Scan Healing	Transfer Abt.	Multi Scan Body	Link Scan Body
							
		TSNSB*	TSSB*	TSSHA*	GSTA/TSTA*	TSMSBC	TSPTB*
	Target	Custom Abt	Transfer Abt.	Link Abt (A)	Link Abt (B)	Multi Ti-Base	Multi/ Multi Ang.
							
		TSPM*/*PM*	GSTA/TSTA*	TSPTB*	TSPL* / TSL*	TSMTB*	GS/TSMA* GS/TS*MA*

* TSPL : Link Abutment with EbonyGold screw
 * TSL : Link Abutment with Ti-screw

List of Digital Library : Exocad

System	Name of library	Scan Source		Target		
		Product Name	Code	Product Name	Code	
SS	Osstem_SS_SBO_CU	New Scan Body	SSNSB*	Custom Abt		SSPM*
	Osstem_SS_SBO_CO			Stock	ComOcta Abt	SSCA*
	Osstem_SS_SBO_COP				ComOcta Plus Abt	SSCAP*
	Osstem_SS_SBO_LI			Link	Link Abt	SSPL*
	Osstem_SS_SB_CU_P	Scan Body	SSSB*	Custom Abt		SSPM*
	Osstem_SS_SB_CO			Stock	ComOcta Abt	SSCA*
	Osstem_SS_SB_COP				ComOcta Plus Abt	SSCAP*
	Osstem_SS_SB_LI			Link	Link Abt	SSPL*
	Osstem_SS_SH_CU	Scan Healing	SSSHA*	Custom Abt		SSPM*
	Osstem_SS_CO_CO	ComOcta Abt	SSCA*	Stock	ComOcta Abt	SSCA*
	Osstem_SS_COP_COP	ComOcta Plus Abt	SSCAP*		ComOcta Plus Abt	SSCAP*

List of Digital Library : Exocad

SS	Scan Source	New Scan body	Scan Body	Scan Healing	ComOcta Abt.	ComOta Plus Abt.	
							
		SSNSB*	SSSB*	SSSHA*	SSCA*	SSCAP*	
SS	Target	Custom Abt	ComOcta Abt.	ComOta Plus Abt.	Link Abt.		
							
		SSPM*	SSCA*	SSCAP*	SSPL*		

List of Digital Library : Exocad



X exocad

System	Name of library	Scan Source		Target		
		Product Name	Code	Product Name		Code
US	Osstem_US_SBO_CU	New Scan Body	USNSB*	Custom Abt		USPM*
	Osstem_US_SBO_CE			Stock	Cement Abt	CAM/CAR/CAW*
	Osstem_US_SBO_LI			Link	Link Abt	USPL*
	Osstem_US_SB_CU	Scan Body	USSB*	Custom Abt		USPM*
	Osstem_US_SB_CE			Stock	Cement Abt	CAM/CAR/CAW*
	Osstem_US_SB_LI			Link	Link Abt	USPL*
	Osstem_US_SH_CU	Scan Healing	USSHA*	Custom Abt		USPM*
	Osstem_US_CE_CE	Cement Abt	CAM/CAR/CAW*	Stock	Cement Abt	CAM/CAR/CAW*

List of Digital Library : Exocad



X exocad

US	Scan Source	New Scan body	Scan Body	Scan Healing	Cement Abt.		
							
		USNSB*	USSB*	USSHA*	CAM/CAR/CAW*		
	Target	Custom Abt	Cement Abt.	Link Abt.			
							
		USPM*	CAM/CAR/CAW*	USPL*			