






# Crestal Approach Sinus-Kit CAS-Kit Drill Guide

● : Required ● : Optional

Fixture selection		Twist Drill	CAS-Drill						Hydraulic Lift & Bone Condensing					CAS-Drill	
			F $\phi$ 4.0		F $\phi$ 4.5		F $\phi$ 5.0		Depth Gauge	Hydraulic Lifter	Bone Carrier	Bone Condenser	Bone Spreader		
Diameter ( $\phi$ )	Fixture protruding height	$\phi$ 2.0	$\phi$ 2.8	$\phi$ 3.1	$\phi$ 3.3	$\phi$ 3.6	$\phi$ 3.8	$\phi$ 4.1							
F $\phi$ 4.0 	0-3 mm	●	●						●	●	●	●	●	●	$\phi$ 3.3
	3-6 mm	●	●	●					●	●	●	●	●	●	-
F $\phi$ 4.5 	0-3 mm	●	●		●				●	●	●	●	●	●	$\phi$ 3.8
	3-6 mm	●	●	●		●			●	●	●	●	●	●	-
F $\phi$ 5.0 	0-3 mm	●	●		●		●		●	●	●	●	●	●	$\phi$ 4.1
	3-6 mm	●	●	●		●		●	●	●	●	●	●	●	-

